

DCDiag

Create a folder such as c:\logfiles\ on the server. Open Windows Explorer and select c:\local drive and create folder called \logfiles\

Open a command prompt command line window and change the working directory to C:\logfiles\ by typing

```
cd\logfiles\
```

and (enter)

Type:

```
dcdiag.exe /f:c:\logfiles\dcdiaglog1.log
```

This will run the DcDiag Tool and save the logfile to the log filename specified.

From:

<https://installconfig.com/> - **Install Config Wiki**

Permanent link:

<https://installconfig.com/doku.php?id=wiki:dcdiag&rev=1509336135>

Last update: **2017/10/30 04:02**

